

Plate Care

Hydro Image

Plate Cleaner and Preserver

Part of CLC SYSTEM

Hydro Image Plate Cleaner contains a proprietary blend of acids which completely desensitize the non-image area of the plate while helping to remove calcium deposits. The pure gum used, along with special humectants and water-attractive film formers, ensure a highly water receptive non-image area.

- Unique blend of acids ensures a desensitized, water receptive
- Non-image area
- Cleans and removes scratches and oxidation.
- Cuts and removes ink quickly
- Removes calcium deposits
- Provides short-term plate storage

Directions of Use:

Pour a liberal amount of Hydro Image Plate Cleaner onto a clean shop towel. Wipe over the entire plate, paying particular attention to the sensitive areas. Wash off with water and run.

Ask your distributor or Allied about other dedicated CLC System Products:

A coordinated pressroom chemistry program designed to control calcium and surfactant contamination across the fountain loop, plates, blankets, rollers, wash system, and fountain solution.



***Important Note: Shake Well Before Use!**

Specifications	
VOC Content lb / gallon	4.5
VOC Content grams / liter	534
Flashpoint	107 °F
Color	White
CA Compliant	No

Available In		
Product #	Size	Case Quantity
10057C	Quart	6